




Modifications for Unique Feet

Everyone has unique feet and individual needs! With this in mind, we've designed Correct Toes to be customizable. Use the chart below to modify your Correct Toes and optimize your toe spacing experience.

Problem	Modification(s)
Skin irritation or rubbing	<ol style="list-style-type: none"> 1. Wear toe socks. 2. Ensure footwear is wide enough. 3. Use one of the following modifications to trim away bothersome material.
Too much spread on pinky toe	<ol style="list-style-type: none"> 1. Cut off half of the pinky spacer. 2. Cut off the full pinky spacer. Shown in GREEN
Correct Toes slide off foot	<ol style="list-style-type: none"> 1. Try toe socks, or a conventional sock over CT (stretch socks' toe seams first). 2. Trim big toe and/or pinky toe pillar, to allow CT to slide further on. Shown in BLUE
Too tight on toe(s)	<ol style="list-style-type: none"> 1. Try a larger size of CT. 2. Cut a small hole on the top and/or bottom. Shown in PURPLE
Bunion Tailor's Bunion	<p>Add a shim in the appropriate CT slot. See YELLOW arrow.</p> <p><i>* Only add shims after 3-4 months of regular Correct Toes use.</i></p> <p>• What's a shim? Take an old shoe liner and cut a piece to fit the slot in your Correct Toes.</p>
Hammertoe	<ol style="list-style-type: none"> 1. Cut a small hole in CT over your hammertoe. Shown in PURPLE 2. Trim away the interspace. Shown in MAGENTA <p><i>* Correct Toes are only effective for hammertoes that are flexible (can be manually straightened).</i></p>



← **Bunion modifications shown in **YELLOW****



Visit www.CorrectToes.com for our Modifications video and detailed information.

→ **Too much spread on pinky toe modification shown in **GREEN****

